

Processing of Autologous Conditioned Serum (ACS) with Orthokine[®]vet irap 60



1

Aseptically withdraw blood using the Orthokine vet irap 60 syringe. Fill syringe, past the soft stop, until plunger locks into place and is flush with the plate and clips. At this point the score mark on the plunger is revealed.

- Apply supplied green luer lock cap and seal securely.
- Ensure the blood and beads are mixed well by 'rocking' the syringe back and forth several times.
- Incubate the Orthokine vet irap 60 syringe at **98.6°F (37°C)** for a **maximum of 24 hours**.



2

After incubation, snap off plunger and centrifuge Orthokine vet irap 60 for 10 minutes at 3700rpm.

- **Ensure the luer lock cap is tightly closed**
- **Load centrifuge symmetrically**



3

After centrifugation, break off the diagonal clips on the guide plate of Orthokine vet irap 60 syringe.

- **Disinfect syringe and port thoroughly**
- **Work on disinfected surface**



4

Take off guide plate by rotating it 90°. Place triangular port on the side of the syringe where serum is the deepest by rotating the base of the syringe plunger. Work carefully to not disturb the clot.



5

Connect cannula to a 35mL syringe. Penetrate the black stopper through accessory port and carefully aspirate ACS, being careful to not aspirate any RBCs.



6

Aliquot ACS through filter/adaptor into 5mL syringes. Close 5mL syringes with a sterile cap. Label each syringe with patient information and date.

- **Store syringes at -18°C (0°F) or colder for up to 7 months**
- **Thaw ACS at room temperature before injecting**

For Technical Support or information, contact:

Dechra Veterinary Products 855-332-9334 or technical.ca@dechra.com

 **Orthokine[®]vet
irap 60**